It is time to Automate the Preparation of Your IFA Slides ...

And Walk Away from the DARK ROOM !

VisiQuantTM

Hyperion ... VisiQuantTM the FIRST TOTALLY AUTOMATED IFA Testing



VisiQuant[™] the FIRST and only AUTOMATED IFA TESTING

- Walk away sample and slide processing
- Titers included with initial reading
- Computerized image processing and storage
- Image comparisons to known ANA patterns
- Customized reports including pattern images
- Archiving of pattern images and results
- LIMS compatible
- Proprietary reagents for optimal performance

Hyperion

14100 SW 136th Street Miami, Florida 33186

> Tel: 305.238.3020 Fax: 305.232.7375

E-mail: hyperion@hyperionclinical.com

www.hyperionclinical.com

VisiQuant™ ANA Test Kit

- ♦ HEp–2 tissue culture cells
- Anti-human IgG, conjugated with unique patented fluorophore
 - La Jolla Blue
 - Near Infrared Excitation and Reading
 - Low auto fluorescence
 - High photostability
- Titers obtained from single initial reading
- FDA cleared for in-vitro diagnostic use
- Additional test kits available soon

